



SPECIFICATIONS

Product Description: **FLEXBAK®**
 Part Number: **7115 SERIES**
 Style: **BACK SUPPORT**

DESCRIPTION:

Flexbak helps reduce back stress and strain with a design so unique, it's patented. Its three-part support is key to a system that provides maximum back support with ease and movement. Because no two bodies are exactly the same, this 9"(22.8 cm) medical-style truss uses two internal support panels above and below the waist surrounded by an elastic third panel. These combine to produce maximum bracing without pinching or riding up and stays deliver firm, continuous lower back and abdominal support. Hand washable.



MATERIAL:

- Spandex style
- Breathable nylon (85.5%)
- Knit material (14.5% Lycra)
- Black

BACK:


- Nine inch wide
- Medical grade truss elastic

CLOSURE:

- Two part hook and loop closures

SUSPENDERS:

- 1"(2.5 cm) nylon removable
- Breakaway
- Fully adjustable
- Non-elastic
- Hook and loop closures



STAYS:

- 2 each molded Polypropylene

TENSION:

- 2 single system pull straps

AVAILABLE:

Part Number	Size	Measurements In (cm)
7115-01	Small	26" (66.4 cm) to 36" (91.4 cm)
7115-02	Medium	36" (91.4 cm) to 48" (121.9 cm)
7115-03	Large	48" (121.9 cm) to 58" (147.3 cm)
7115-04	X-Large	58" (147.3 cm) to 68" (172.7 cm)